

Polarity Index	Solvent	BP (°C)
0.1	heptane	98.43
0.2	cyclohexane	80.72
2.4	toluene	110.62
2.5	t-butyl methyl ether	55.2
3.1	dichloromethane	39.75
3.5	ethylene dichloride	83.48
3.9	<i>n</i> -butyl alcohol	117.5
4.0	tetrahydrofuran	66
4.1	chloroform	61.15
4.4	ethyl acetate	77.11
4.7	methyl ethyl ketone (MEK)	79.64
4.8	dioxane, 1,4-	101.32
5.1	acetone	56.29
5.1	methanol	64.7
5.3	pyridine	115.25
5.8	acetonitrile	81.6
7.2	dimethyl sulfoxide	189
10.2	water	100

References: Journal of Chromatography A, Volume 92, Issue 2, 22 May 1974, Pages 223-230.

[illegible]

■ Water miscible ■ Miscible ■ Immiscible